



Impression Trays

Made of 18/8 Stainless Steel, Solid, Set of 10



All The Above Impression Trays are Available in These 5 Sizes



UPPER		U1	U2	U3	U4	U5
	A	62	60	53	54	52
	B	83	75	73	69	68


All Dimension are in mm



LOWER		L1	L2	L3	L4	L5
	A	63	61	57	60	53
	B	82	81	78	74	71

All Dimension are in mm

Perforated, Set of 10



Impression Trays

Made of 18/8 Stainless Steel, Solid, Full Denture Set of 8




All The Above Impression Trays are Available in These 4 Sizes

UPPER		XL	L	M	S
	A	70	67	64	60
LOWER		XL	L	M	S
	A	64	64	60	57
UPPER	B	79	76	70	73
	B	86	70	73	70

All Dimension are in mm


All Dimension are in mm



Perforated Full Denture Set of 8

Impression Trays

Made of 18/8 Stainless Steel, Solid, Set of 6




All The Above Impression Trays are Available in These 4 Sizes

UPPER		U1	U2	U3
	A	67	64	57
	B	81	71	70

All Dimension are in mm

LOWER		L1	L2	L3
	A	60	57	54
	B	81	79	76

All Dimension are in mm




Perforated Set of 6

Impression Trays

Made of 18/8 Stainless Steel, English Pattern Set of 8




All The Above Impression Trays are Available in These 4 Sizes



UPPER	Size	A	B	C
	0	67	50	12
	1	72	51	13
	2	78	55	15
	3	80	58	15

All Dimension are in mm




LOWER	Size	A	B	C
	0	70	50	14
	1	77	51	15
	2	79	55	15
	3	80	60	17


All Dimension are in mm

Impression Trays

Made of 18/8 Stainless Steel, Solid, Full Denture Set of 8



All The Above Impression Trays are Available in These 4 Sizes




UPPER		XL	L	M	S
	A	70	67	64	60
	B	79	76	70	73

All Dimension are in mm



LOWER		XL	L	M	S
	A	64	64	60	57
	B	86	70	73	70


All Dimension are in mm



Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8

Impression Trays

Made of 18/8 Stainless Steel, Solid, Full Denture Set of 8




Upper Sizes: 0,1,2,3

Depth of Tray: 15,9 mm, Height of Plate: 12,5 mm

Lower Sizes: 0,1,2,3



Depth of Tray: 12,9 mm



Made of 18/8 Stainless Steel, Perforated Full Denture Set of 8

Impression Trays

Impression Trays For Crown & Bridge Work



2931-A
For upper right jaw
For lower left jaw

2931-B
For upper left jaw
For lower right jaw

Available in 4 sizes

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm

2931-C
For upper right jaw
For lower left jaw

2931-D
For upper left jaw
For lower right jaw

Available in 4 sizes

- 1.....44 mm
- 2.....55 mm
- 3.....60 mm
- 4.....65 mm